

TECNO 3 "CCC"

COMPLETE CLOSED CIRCUIT

New

BIOS.C

Line



OZONOTHERAPY

Closed-loop medical device intended for ozonated self-motor therapy that consists of blood sampling, followed by treatment with gaseous oxygen/ozone mixture and subsequent reinfusion to the patient.

FIELDS OF APPLICATION:

Treatment of disease states in which stimulation of the immune system and/or increased tissue oxygenation is required.



TECHNICAL SPECIFICATIONS

PHthalates FREE sterile single-use closed-loop medical device with pre-connected withdrawal and re-infusion line, drip chamber with 200 microns filter and with needle-free connections for needle-free use.

ADVANTAGES:

- No contact with environment outside the device
- Emergency by-pass system in case of re-infusion line failure
- Proximal connection to the patient for additional drug withdrawals or infusions
- Connection for insufflation of Ozone-Oxygen mixture equipped with one-way valve to prevent any spillage and anti-bacterial filter to protect against external contamination.



OZONE THERAPY KITS AVAILABLE

COD.	DESCRIPTION	Q.TÀ/CONF.
TEC865	TECNO3 CCC with built-in variable dispenser (ACDA anticoagulant)	15 pcs
TEC876	TECNO3 CCC with prefilled bag with 20ml of ACDA anticoagulant	20 pcs
TEC876/35	TECNO 3CCC with prefilled bag with 35ml of ACDA anticoagulant	20 pcs


NEW BIOS.C LINE

CCP876	TECNO3 CCC with prefilled bag with 20ml of ACDA anticoagulant	<u>BIOS.C</u>	20 pcs
CCP876/35	TECNO 3CCC with prefilled bag with 35ml of ACDA anticoagulant	<u>BIOS.C</u>	20 pcs

CE 0051

BIOS.C LINE

The **TECNO3 "CCC" BIOS.C** device is made with a coating treatment of surfaces in contact with blood with nonthrombogenic characteristics consisting of a polymer film with prosthetic groups of phosphorylcholine; thus, a bio-mimetic structure simulating the structure of cell surfaces is obtained. This mechanism of biomimicry makes it possible to reduce the adverse reactions that are established when blood vessels are in contact with poorly biocompatible foreign materials.



ACCESSORIES

EO51P	Sampling scale	1 pc
-------	----------------	------



Tecnoline S.p.A.
 Registered Office Via Gelatti 11/A - 41033 Concordia sulla Secchia - Modena - Italy
 Headquarters Via Gelatti 11/A and 18 - 41033 Concordia sulla Secchia - Modena - Italy
 T 0535 40312 M info@tecnolinemed.com