

IGIENITEC

Complete closed device



IGIENITEC

It is a device that can prevent environmental and skin contamination by fecal "soiling."

PARTICULARLY SUITABLE:

- HOSPITAL USE:
 - Postoperative management to ensure maximum cleanliness;
- HOME USE:
 - Bed-restrained subjects or those with poor physical and mental self-sufficiency and urinary and fecal incontinence.
 - Daily fecal continence control in active and self-sufficient people but with diarrheal bowel disease and insufficient sphincter control.

IGIENITEC consists of a flexible wafer that has an opening in the center. One side of the wafer adheres to the skin around the anus, and the other side is connected to a collection pouch. The bag has a closable attachment at the lower end through which a catheter can be inserted to perform enteroclysms, with maximum hygienic precaution and environmental isolation.

Access is closed by a ventilated capsule that allows flatulence to escape by filtering out foul odor and preventing distension of the pouch by gas entrapment.


Together with the adhesive surface that taps into the anal rhyme is made of totally latex-free material, it unlike existing devices IS ABSOLUTELY NOT INVASIVE and does not require the insertion of a balloon catheter.

Patented Product: EP 1 768 624 B1

AVAILABLE DEVICES

COD.	Description	Q.TÀ/CONF.
TEC353608BL	IGIENITEC	25 pcs



 Tecnoline S.p.A.
Registered Office Via Gelatti 11/A - 41033 Concordia sulla Secchia - Modena - Italy
Headquarters Via Gelatti 11/A and 18 - 41033 Concordia sulla Secchia - Modena - Italy
T (+39) 0535 40312 M info@tecnolinemed.com

