

SCAN ME



## OZONE THERAPY

Cod. TEC773

Medical device with a closed circuit, intended for both Major Ozonated Autohemotherapy (GAEI)—which involves blood collection, followed by treatment with an oxygen/ozone gas mixture and subsequent reinfusion—and for the ozonation of saline solution.

## TECHNICAL SPECIFICATIONS

### TEC876 TEC876/35

Sterile closed-circuit medical device  
PHTHALATES FREE, with pre-connected blood collection and reinfusion lines, drip chamber with 200µ blood filter, needle-free connections for needleless use, and graduated bag pre-filled with ACD-A anticoagulant. All components in contact with ozone, particularly the ozone inlet, are made of ozone-resistant medical-grade PVC, phthalate-free and latex-free, to ensure maximum chemical compatibility and safety during treatment.

STERILE R

### ADVANTAGES AND COMPOSITION TEC876

- ▶ Dual pathway for ozone inlet and outlet
- ▶ Safety valves for ozone inflow and outflow
- ▶ Emergency bypass in case of reinfusion line interruption
- ▶ Device pre-filled with ACD-A

## TECHNICAL SPECIFICATIONS

### TEC773

TEC773 is a closed-circuit medical device  
PHTHALATES FREE, designed to ensure safety, convenience, and reliability in the systemic administration of ozone.  
Ideal for clinical and outpatient use, it offers a complete and ready-to-use set for efficient therapy management. All components in contact with ozone, particularly the ozone inlet, are made of ozone-resistant medical-grade PVC, phthalate-free and latex-free, to ensure maximum chemical compatibility and safety during treatment.

STERILE EO

### ADVANTAGES AND COMPOSITION TEC773

- ▶ Dual pathway for ozone inlet and outlet
- ▶ Safety valves for ozone inflow and outflow
- ▶ Emergency bypass in case of reinfusion line interruption
- ▶ 2 extensions, one for ozone inlet and one for outlet, to connect to the generator
- ▶ Universal spike for device filling

## ORDERING GUIDE

| CODE      | DESCRIPTION                                   | QUANTITY |
|-----------|---|----------|
| TEC876    | TECNO3 pre-filled with 20ml of ACD-A          | 20 PZ    |
| TEC876/35 | TECNO3 pre-filled with 35ml of ACD-A          | 20 PZ    |
| TEC773    | TECNO3 Advanced system for systemic ozonation | 20 PZ    |



Tecnoline S.p.A.  
Registered Office Via Gelatti 11/A - 41033 Concordia sulla Secchia - Modena - Italy  
Headquarters Via Gelatti 11/A and 18 - 41033 Concordia sulla Secchia - Modena - Italy  
T (+39) 0535 40312 M info@tecnolinemed.com



Rev.5 17/07/2024