

SCAN ME



OZONE THERAPY

Cod. TEC876 - TEC876/35

Medical device with a closed circuit, intended for Major Ozonated Autohemotherapy (GAEI), which involves blood collection, followed by treatment with an oxygen/ozone gas mixture, and immediate reinfusion.

TECHNICAL SPECIFICATIONS

TEC876 TEC876/35

Sterile closed-circuit medical device
PHTHALATES FREE, with pre-connected blood collection and reinfusion lines, drip chamber with 200µ blood filter, needle-free connections for needleless use, and graduated bag pre-filled with ACD-A anticoagulant. All components in contact with ozone, particularly the ozone inlet, are made of ozone-resistant medical-grade PVC, phthalate-free and latex-free, to ensure maximum chemical compatibility and safety during treatment.

STERILE R

ADVANTAGES AND COMPOSITION TEC876

- ▶ Dual pathway for ozone inlet and outlet
- ▶ Safety valves for ozone inflow and outflow
- ▶ Emergency bypass in case of reinfusion line interruption
- ▶ Device pre-filled with ACD-A

TECHNICAL SPECIFICATIONS

TEC773

TEC773 is a closed-circuit medical device
PHTHALATES FREE, designed to ensure safety, convenience, and reliability in the systemic administration of ozone.
Ideal for clinical and outpatient use, it offers a complete and ready-to-use set for efficient therapy management. All components in contact with ozone, particularly the ozone inlet, are made of ozone-resistant medical-grade PVC, phthalate-free and latex-free, to ensure maximum chemical compatibility and safety during treatment.

STERILE EO

ADVANTAGES AND COMPOSITION TEC773

- ▶ Dual pathway for ozone inlet and outlet
- ▶ Safety valves for ozone inflow and outflow
- ▶ Emergency bypass in case of reinfusion line interruption
- ▶ 2 extensions, one for ozone inlet and one for outlet, to connect to the generator
- ▶ Universal spike for device filling

ORDERING GUIDE

CODE	DESCRIPTION	QUANTITY
TEC876	TECNO3 pre-filled with 20ml of ACD-A	20 PZ
TEC876/35	TECNO3 pre-filled with 35ml of ACD-A	20 PZ
TEC773	TECNO3 Advanced system for systemic ozonation	20 PZ



Tecnoline S.p.A.
Registered Office Via Gelatti 11/A - 41033 Concordia sulla Secchia - Modena - Italy
Headquarters Via Gelatti 11/A and 18 - 41033 Concordia sulla Secchia - Modena - Italy
T (+39) 0535 40312 M info@tecnolinemed.com



Rev.5 17/07/2024